

Basic Drapes


Dina-Hitex drapes are available in various sizes, materials and apertures.

Simple Drapes

May be used as table covers or wrapping sheets.

Size (cm)	Ref SMS	Ref DS2	Ref DS3	Ref 2SBL	Ref FS2	Ref FS3	Ref PF2
34 x 34		20-0169-S					
50 x 50	10-0101-S	20-0101-S	30-0101-S	51701-S	11-0101-S	22-0101-S	25-0101-S
50 x 70	10-0102-S	20-0102-S	30-0102-S	51702-S	11-0102-S	22-0102-S	25-0102-S
75 x 75	10-0103-S	20-0103-S	30-0103-S	51703-S	11-0103-S	22-0103-S	25-0103-S
75 x 90	10-0104-S	20-0104-S	30-0104-S	51704-S	11-0104-S	22-0104-S	25-0104-S
100 x 120	10-0132-S	20-0132-S	30-0132-S	51717-S	11-0132-S	22-0132-S	25-0232-S
150 x 180	10-0106-S	20-0106-S	30-0106-S	51783-S	11-0106-S	22-0106-S	25-0106-S
200 x 240	10-0133-S	20-0133-S	30-0133-S	51721-S	11-0133-S	22-0133-S	*

* PF2 max possible size 150 x 240 cm




Surgical Drapes with AT on the Edge

Maximum length of AT is defined by the length of the drape edge.

Size (cm)	AT (cm)	Ref SMS	Ref DS2	Ref DS3	Ref 2SBL	Ref FS2	Ref FS3	Ref PF2
50 x 50	50	10-0201-S	20-0201-S	30-0201-S	51772-S	11-0201-S	22-0201-S	25-0201-S
50 x 70	70	10-0202-S	20-0202-S	30-0202-S	51752-S	11-0202-S	22-0202-S	25-0202-S
75 x 75	75	10-0203-S	20-0203-S	30-0203-S	51742-S	11-0203-S	22-0203-S	25-0203-S
75 x 90	90	10-0204-S	20-0204-S	30-0204-S	51735-S	11-0204-S	22-0204-S	25-0204-S
100 x 120	120	10-0232-S	20-0232-S	30-0232-S	51740-S	11-0232-S	22-0232-S	25-0232-S
150 x 180	80	10-0206-S	20-0206-S	30-0206-S	51812-S	11-0206-S	22-0206-S	25-0206-S
200 x 240	80	10-0233-S	20-0233-S	30-0233-S	51738-S	11-0233-S	22-0233-S	*


* PF2 max size 150 x 240 cm



Surgical Drapes with Aperture

Apertures may be usually of either round, oval, square, rectangular, diamond or pentagon shape.

Size (cm)	Aperture Size (cm)	Aperture Shape	Ref DS2	Ref DS3	Ref 2SBL	Ref FS2	Ref FS3
50 x 70	7	round	20-0301-S	30-0301-S	15-0311-S	11-0301-S	22-0301-S
50 x 70	10	round	20-0302-S	30-0302-S	15-0312-S	11-0309-S	22-0302-S
75 x 90	7	round	20-0303-S	30-0303-S	15-0305-S	11-0303-S	22-0303-S
75 x 90	10	round	20-0304-S	30-0304-S	15-0317-S	11-0316-S	22-0304-S
75 x 110	10	round	20-0305-S	30-0305-S	15-0315-S	11-0305-S	22-0305-S
50 x 50	6 x 8	oval	20-0308-S	30-0308-S	15-0303-S	11-0308-S	22-0308-S
50 x 70	6 x 8	oval	20-0314-S	30-0314-S	15-0314-S	11-0314-S	22-0314-S
100 x 150	6 x 8	oval	20-0306-S	30-0306-S	15-0316-S	11-0306-S	22-0306-S
160 x 200	12	round	20-0366-S				



BARRIER SYSTEM FOR SURGICAL THEATRES

Sterile Surgical Drapes and Sets


Basic Drapes

Surgical Drapes with Adhesive Aperture

The adhesive area is made by AT placed around the aperture.

Size (cm)	Aperture Size (cm)	Aperture Shape	Ref DS2	Ref DS3	Ref 2SBL	Ref FS2	Ref FS3
50 x 70	7	round	20-0401-S	30-0401-S	15-0403-S	11-0401-S	22-0401-S
50 x 70	10	round	20-0402-S	30-0402-S	15-0404-S	11-0422-S	22-0402-S
75 x 90	7	round	20-0403-S	30-0403-S	15-0419-S	11-0403-S	22-0403-S
75 x 90	10	round	20-0404-S	30-0404-S	15-0409-S	11-0404-S	22-0404-S
75 x 110	10	round	20-0405-S	30-0405-S	15-0420-S	11-0423-S	22-0405-S
90 x 120	7	round	20-0422-S	30-0406-S	15-0421-S	11-0406-S	22-0406-S
120 x 150	10	round	20-0407-S	30-0407-S	15-0422-S	11-0407-S	22-0407-S
45 x 75	7	round	20-0409-S	30-0409-S	15-0423-S	11-0409-S	22-0409-S
150 x 200	7	round	20-0446-S	30-0446-S	15-0424-S	11-0424-S	22-0424-S
34 x 34	10	round	20-0512-S				
2 x (34 x 34) (2 part drape)	2 x (6 x 8)	2 x U	20-0513-S				


Aperture options



Cannular and Catheterisation Drapes

Drapes with aperture and radial perforation made of either laminated DS2 or transparent PE.

Size (cm)	Aperture Size (cm)	Accessories	Ref DS2	Ref PE	Item
75 x 100	10	AT	501-005-S	-	Cannular drape
75 x 90	10	-	701-181-S	701-185-S	Drape with aperture
90 x 120	10	-	701-188-S	701-192-S	
120 x 150	10	-	701-195-S	-	
75 x 90	10	AT	701-182-S	701-186-S	Drape with adhesive aperture
90 x 120	10	AT	701-189-S	701-193-S	
120 x 150	10	AT	701-196-S	-	
75 x 90	10	PE+AT	701-183-S	-	Drape with PE and adhesive aperture
90 x 120	10	PE+AT	701-190-S	-	
120 x 150	10	PE+AT	701-197-S	-	
75 x 90	10	AT+AT	701-184-S	701-187-S	Drape with adhesive aperture and adhesive edge
90 x 120	10	AT+AT	701-191-S	701-194-S	
120 x 150	10	AT+AT	701-198-S	-	



Absorbent Surgical Table Covers

Size (cm)	Material	Reinforcement (cm)	Code	Item
100 x 220	FS2	SM-M	50891	Super absorbent table cover
140 x 220	FS2	SM-M	50892	Super absorbent table cover
100 x 240	SMS35g	FSM80	50894	Absorbent table cover
140 x 240	SMS35g	FSM80	50895	Absorbent table cover

