


Cardiovascular Surgery

Angio drapes are usually equipped with high absorbent reinforced area to retain leaking blood and with 2 or 4 apertures with IF or AT. Some drapes have PE transparent foil on one or both sides. Cardiovascular drapes can be equipped with tube holders, fixing tapes and collection pouches.


Angio drape with 2 PE transparent sides and 2 adhesive apertures.

Code	501-241-S
Material	DS2
Size (cm)	210 x 300
Aperture (cm)	2x (6 x 8)
Accessories	AT




Angio drape, reinforced, with 2 PE transparent sides and 2 apertures in IF.

Code	501-239-S
Material	DS2
Size (cm)	220 x 365
Reinforcement (cm)	SM 65 x 77
Aperture (cm)	2x (12)
Accessories	IF with aperture 8 cm



Angio drape, reinforced, with 2 PE transparent sides and 2 apertures with IF.

Code	501-240-S
Material	SMS
Size (cm)	220 x 330
Reinforcement (cm)	FSM80 75 x 120
Aperture (cm)	2x (12)
Accessories	IF




BARRIER SYSTEM FOR SURGICAL THEATRES

Sterile Surgical Drapes and Sets

Cardiovascular Surgery


Angio drape, reinforced, with 2 PE transparent sides and 4 apertures in IF.

Code	501-242-S
Material	SMS+PE
Size (cm)	240 x 345
Reinforcement (cm)	FSM 122 x 130
Aperture (cm)	2x (14), 2x (9 x 12)
Accessories	IF with aperture 2x (11), 2x (7 x 9)



Cardiovascular drape with IF, 2 pouches and 3 tube holders.

Code	503-100-S
Material	DS3
Size (cm)	280 x 340
Aperture (cm)	32 x 40
Accessories	IF, holders, 2 pouches DS2



Pectoral drape, reinforced, with IF, a holder and 2 pouches.

Code	503-113-S
Material	SMS
Size (cm)	260 x 360
Reinforcement (cm)	FSM80 100 x 200
Aperture (cm)	32 x 40
Accessories	IF, holder, 2 pouches FSM80

