

BARRIER SYSTEM FOR SURGICAL THEATRES

Sterile Surgical Drapes and Sets

Delivery and Newborn Sets

Sets are applied on women during labour, suitable for newborns as well. Using disposable products reduces risk of infection.

Common items used in these sets:

Bed sheet	SB, 2SBL
Drapes	variety of types and materials
Impervious protective apron	PE
Impervious sleeve cover	PE
Isothermic foil	
Limb covers	SMS
Table cover	FS2, 2SBL
Protective gown	SB
Surgical gown PRIMA	SMS
Instruments:	Scalpel, scissors sharp, umbilical clamp
Umbilical tie	Plastics
Abdominal sponge	Gauze, NW
Absorbent pad	ONTEX
Absorbent towel	Vis, tissue
Compress	Gauze, NW
Cotton sheet	
Sanitary towel	
Twisted tampon	Gauze
Gloves	Latex



Delivery and New-born Sets

Delivery Set 60550-44-S

Component	Material	Size (cm)	Pcs
Wrapping sheet	FS2	75 x 80	1
Instruments table drape	2SBL	80 x 100	1
Absorbent towel	VIS40	20 x 30	6
Umbilical tie		1 x 15	1
Umbilical clamp			2
Absorbent towel	Spunlace	70 x 80	1
Absorbent pad Santex	Cellulose	60 x 90	1
PE-Apron, white	PE	80 x 125	1
Surgical drape	2SBL	75 x 130	1

Delivery Set 60550-52-S

Component	Material	Size (cm)	Pcs
Instruments table drape	2SBL	80 x 100	1
Surgical drape	2SBL	75 x 130	1
Surgical gloves, powdered, 1 pair	Latex	Size 8	2
PE-Apron, white	PE	80 x 125	2
Scalpel type 11			1
Absorbent towel	VIS40	20 x 30	2
Umbilical tie		1 x 15	3
Umbilical clamp			2
Sanitary towel			2
Wrapping sheet	FS2	60 x 60	1
Surgical drape	2SBL	60 x 60	1
Isothermic foil		70 x 90	1
Cotton sheet		70 x 70	2
Suction catheter, 6ch			1

Delivery Set 60550-57-S

Component	Material	Size (cm)	Pcs
Wrapping sheet	FS2	75 x 80	1
Absorbent pad Santex	Cellulose	60 x 90	1
Umbilical clamp			2
Scissors-sharp	Stainless	11,5	1
Examination gloves, powdered	Latex	M	2
Examination gloves, powdered	Latex	L	2
Umbilical tie		1 x 15	2
Abdominal drape	NW	60 x 60	1
Gauze compress	Gauze	10 x 10	10
Tampon twisted	Gauze	30 x 30	10