

DISPOSABLE CLOTHING FOR HOSPITALS

Hygienic Sheets, Bed Sheets and Bedding Sets

Hygienic Sheets Hitex

Made of non-woven PP hydrophobic SB, non-irritant and non-allergenic material. Used as a protection of examination chairs and beds during outpatient treatment. Available in rolls from 45 up to 100 cm wide and either 100 m or 300 m long. Colour white and blue.


Width (cm)	Ref white 23 g 100 m	Ref white 23 g 300 m	Ref blue 28 g 100 m	Ref blue 28 g 300 m	Ref blue 17g 100 m	Ref blue 17g 300 m	Ref white 17g 100 m	Ref white 17g 300 m
45	50100	50400						
50	50101	50401	50121	50421			50101-B	50401-B
55	50102	50402	50122	50422	50102-M	50402-M		
60	50103	50403	50123	50423	50103-M	50403-M	50103-B	50403-B
70	50104	50404						
80	50105	50405						
100	50106	50406						



Hygienic Sheets Hitex Pro

Made of highly absorbent 2ply laminated material, consists of cellulose and water proof PE film, provides maximum protection and comfort. Available in rolls, cross perforations ensure an easy tearing of the sheets.

Ref	HitexPro Colour	Sheet size (cm)	Sheets per roll
58102	Green	33 x 50	50
58103	Blue	33 x 50	50
58212	Green	51 x 50	80
58213	Blue	51 x 50	80
58222	Green	51 x 80	50
58223	Blue	51 x 80	50
58232	Green	51 x 160	25
58233	Blue	51 x 160	25



Roll Towel

Made of absorbent MB. Available in rolls 108 m long, cross perforations ensure an easy tearing of the sheets.

Ref	Sheet size (cm)	Sheets per roll
96001	32 x 32	338



Towel

Made of airlaid cellulose.

Ref	Size (cm)
P-3800-001	30 x 33cm



Hygienic Sheets, Bed Sheets and Bedding Sets

Bed Sheets

These bed sheets are cost-effective substitution of linen or cotton bed sheets during outpatient treatment or short term hospitalization. Made of non-irritant non woven textile either SB or SMS or water proof laminated material 2SBL and PE. 2SBL and PE provide an effective barrier against fluids penetration.

Ref	Size (cm)	Material	Colour
50801	80 x 220	SB23	White
50802	100 x 220	SB23	White
50807	160 x 240	SB40	White
50811	80 x 220	SB23	Green
50812	100 x 220	SB23	Green
50817	160 x 240	SB40	Green
50854	100 x 220	PE + elastic band	White
50801-SMS	80 x 220	SMMMMS 50	Blue
50802-SMS	100 x 220	SMMMMS 50	Blue
50801-SMS-35	80 x 220	SMS 35	Blue
51705	80 x 100	2SBL	Blue
51708	100 x 150	2SBL	Blue
51709	100 x 200	2SBL	Blue
51712	120 x 150	2SBL	Blue
51715	150 x 200	2SBL	Blue

Stretcher's Pad Cover

Hygienic ambulance stretcher pad cover made of white SB. Fitted with open pleats for easy fastening on either side of the stretcher.

Ref	Material	Size (cm)
50501DH	SB23	60 x 190
50502DH	SB40	60 x 190

Bedding Set Economic

This Set is cost-effective substitution of linen or cotton bed clothes during outpatient treatment or short term hospitalization. Made of non-irritant SB material.

Item	Ref	Colour
Bedding set Economic: consists of pillow case 70 x 80 cm, quilt case 155 x 200 cm, bed sheet 160 x 240 cm	59001	White
Bedding set Economic: consists of pillow case 70 x 80 cm, quilt case 155 x 200 cm, bed sheet 160 x 240 cm	59007	Green
• Quilt case 155 x 200 cm	56201	White
• Quilt case 155 x 200 cm	56202	Green
• Pillowcase 70 x 80 cm	55201	White
• Bed sheet 160 x 240 cm	50807	White