


Surgical Covers and Accessories

Table Covers

Made of impermeable material FS2 or 2SBL.


Size (cm)	Ref FS2	Ref 2SBL
50 x 70	11-0102-S	51702-S
80 x 100	11-0136-S	51705-S
80 x 120	11-0137-S	51706-S
100 x 120	11-0132-S	51717-S
100 x 150	11-0116-S	51708-S
100 x 200	11-0174-S	51709-S
150 x 200	11-0117-S	51715-S



MAYO Covers

Made of impermeable blue PE material, non-transparent, fitted with anti-slip surface and absorbent reinforcement. The size of the sack shaped cover is 80 x 145 cm.

Reinforced area (cm)	Reinf. Material	Ref
50 x 80	SB28	40-0101-S
70 x 80	SB28	40-0105-S
80 x 145	SB28	40-2001-S
50 x 145	DS2	40-0201-S
70 x 80	DS2	40-0205-S
80 x 145	DS2	40-2005-S
50 x 145	DS2	40-0209-S telescopic folding



Limb Covers

Sack shaped and usually made of laminated absorbent viscose – DS2. Used for easy limb covering during orthopaedic surgeries.

Size (cm)	Ref
30 x 60	50-0201-S
30 x 80	50-0202-S
40 x 120	50-0203-S



URO and GYN Covers

Usually made of hydrophobic SMS material. Used for easy leg covering during transurethral or transvaginal surgeries.

Size (cm)	Ref
80 x 100	50-0102-S
80 x 120	50-0104-S



BARRIER SYSTEM FOR SURGICAL THEATRES


Set's Accessories

Surgical Covers and Accessories

Dinafix

Adhesive tube and cable holder on Velcro base.

Size (cm)	Ref
2,5 x 10	78-0117-S
2,5 x 15	78-0105-S
2,5 x 20	78-0101-S
2,5 x 25	78-0110-S
2,5 x 30	78-0106-S



AT - Adhesive Tape

Adhesive tape with absorbent viscose on the top, used for tube fastening.

Size (cm)	Ref
5 x 30	70-0113-S
10 x 50	70-0101-S
10 x 70	70-0107-S
15 x 20	70-0105-S
15 x 50	70-0103-S



Incise Foils

Adhesive polyurethane sterile foil used for covering operating area, reduces risk of wound contamination.


Size (cm)	Ref
10 x 15	76-0124-S
10 x 20	76-0113-S
15 x 20	76-0101-S
15 x 25	76-0102-S
35 x 40	76-0103-S
20 x 30	76-0104-S
45 x 50	76-0105-S
20 x 45	76-0106-S
25 x 30	76-0107-S
25 x 40	76-0108-S
25 x 50	76-0109-S
40 x 50	76-0110-S
40 x 65	76-0111-S
50 x 80	76-0112-S

Device Covers

Made of PE foil, equipped with an elastic band for easy fastening. Protects device from contamination and keeps the surface sterile.


Round Covers

Size (cm)	Ref
42,5 x 42,5	71-0106-S
50 x 60	71-0125-S
60 x 100	71-0115-S
60 x 120	71-0116-S
65 x 85	71-0103-S
70 x 70	71-0117-S
80 x 120	71-0118-S
80 x 160	71-0119-S
85 x 85	71-0102-S
85 x 95	71-0120-S
95 x 130	71-0121-S
100 x 120	71-0122-S
100 x 125	71-0104-S
105 x 105	71-0108-S
120 x 120	71-0111-S
140 x 140	71-0105-S
180 x 200	71-0107-S



Flat Covers

Size (cm)	Ref
70 x 80	75-0101-S
75 x 75	75-0102-S
80 x 150	75-0103-S
85 x 90	75-0104-S
90 x 100	75-0105-S
90 x 160	75-0106-S
100 x 115	75-0107-S
100 x 160	75-0108-S
120 x 150	75-0109-S



Camera Covers

Ref	Size (cm)	Accessories
80-0101-S	14 x 250	Fixing tape, card cuffs
80-0105-S	14 x 145	Fixing tape, card cuffs
80-0109-S	14 x 250	Elastic tip, card cuffs



80-0101/80-0105



80-0109

BARRIER SYSTEM FOR SURGICAL THEATRES

Set's Accessories

Device Covers

Cord Covers

Ref	Size (cm)	Accessories
80-0103-S	8 x 120	Fixing tape
80-0108-S	8 x 120	Card cuff



80-0103



80-0108

Ultra Sound Covers

Ref	Size (cm)	Accessories
80-0133-S	20 x 120	Adhesive tape, telescopic folding
80-0132-S	20 x 250	Adhesive tape, telescopic folding



Drawstring Cover

2-Drawstring bag is made of the PE film 40 µm.

Ref	Size (cm)	Accessories
80-0135-S	70 x 80	2 drawstring sides



C-Arm Covers - Universal


Made of transparent PE foil.

Designed as a one piece tube cover with fixing tapes or a three-part set.

The three-part set is composed of two 71-0111 covers and a PE sheet with 4 fixing tapes.

Fixing tapes enable easy fastening around C-arm.

Ref	Size (cm)	Accessories
40-0334-S	120 x 180	Fixing tapes
40-0335-S	120 x 260	Fixing tapes
40-0337-S	120 x 200	three-part set




Liquid Collection Pouches

Adhesive Pouches

Designed to collect fluid from the surgery area.
Pouches are made of strong PE transparent foil.
Equipped with adhesive tape at the back side which ensures easy fastening to the drape.
Front part of the pouch is equipped with a forming brace.


1-Compartment Pouch

Size (cm)	Ref
22 x 29	60-0101-S
28 x 44	60-0102-S
28 x 55	60-0103-S




2-Compartment Pouch

Size (cm)	Ref
25 x 28	60-0124-S
28 x 44	60-0125-S



3-Compartment Pouch


Size (cm)	Ref
28 x 44	60-0149-S



Pouch for Neurosurgery

Adhesive tape and semi oval aperture on upper rear side edge enables perfect adjustment to the head shape.
The front part with forming brace enables free fluid access into the pouch.


Size (cm)	Ref
30 x 32	60-0111-S



Pouches with Filter

The rear side of the pouches is equipped with a 5 cm wide adhesive tape.
The front part of the pouch is equipped with a forming brace.
Conical drainage port for tube setting (6 - 10 mm in diameter) is equipped with a filter to prevent from clogging.

Ref	Size (cm)	Accessories
65-0101-S	50 x 50	Drainage, filter
65-0102-S	50 x 60	Drainage, filter
65-0106-S	60 x 90	Drainage, filter



BARRIER SYSTEM FOR SURGICAL THEATRES

Set's Accessories

Liquid Collecting Pouches

Arthroscopic Pouches - Knee, Shoulder

Pouches are designed to collect fluid during knee or shoulder arthroscopy. They are made of strong PE foil, equipped with a conical drainage port for tube setting (6 - 10 mm in diameter) and a forming brace. Tight shape around the limb is ensured by 2 flexible elastic membranes with oval aperture 6 x 8 cm.

Ref	Size (cm)	Accessories
60-0122-S	60 x 100	2 elastic membranes



Caesarean Pouches

Made of strong transparent PE foil. Pouches are equipped with an adhesive tape around the aperture. The front part has either a forming brace or a PUR foam brace to retain fluids safely. Pouches are accessoried with a conical drainage port for tube setting (6 - 10 mm in diam.) on both sides or with a sleeve outlet.

Ref	Size (cm)	Accessories
60-0106-S	80 x 100	PUR foam
60-0107-S	80 x 100	Plastic forming brace
60-0145-S	110 x 110	PE sleeve outlet



Pouches for Gynaecology

They are made of strong PE foil as either 1-compartment or 2-compartment pouches. All pouches are equipped with a forming brace in front and an adhesive tape on rear side to be attached to the drape. Pouches are fitted with a conical drainage port or sleeve outlet and volume indicator.

Ref	Size (cm)	Accessories
60-0146-S	60 x 150	1 compartment, without scale, sleeve outlet
60-0153-S	60 x 100	1 compartment, with scale and drainage port
60-0148-S	60 x 100	2 compartments, with scale and drainage port